


## TiO<sup>2</sup> Surgical Face Mask - Type III



	LightCleanse TiO <sup>2</sup> Surgical Face Mask -Type III	Conventional Surgical Face Mask
Bacterial Filtration Efficiency	Yes	Yes
Particulate Filtration Efficiency	Yes	Yes
Immediate Elimination of Bacteria on Contact	Yes	No
Immediate Elimination of Virus on Contact	Yes	No
Effective even when wet	Yes	No
Prevents Cross Contamination Upon Contact	Yes	No
Easy & Safe Disposal	Yes	No

LightCleanse TiO <sup>2</sup> Surgical Face Mask -Type III
3 Ply
Inner Facing - Layer Polypropylene Thermal Bonded Fabric
Outer Facing - Layer Polypropylene Thermal Bonded Fabric
Filter - Polypropylene Non-Woven
Side Tapes - Polyester Non-Woven
Ear Loop - Polyurethane & Polyester
Nose Clamp - Polyethylene Covered Steel Wire (4.3mm wide and core dia. 0.6mm)

LightCleanse TiO<sup>2</sup> Surgical Face Mask -Type III meets ASTM F1201 & EN14683 Standards